



## Halo Disinfection System®

# The Most Effective and Affordable Disinfection Solution for Emergency Vehicles

GermS found inside the vehicles used to transport injured or ill people can spread infections to the next occupant and EMT professionals. The very real possibility of cross-contamination arises when patients are transferred between healthcare facilities and nursing homes, and can impact the health of hundreds of patients and staff. When it comes to the effectiveness of vehicle disinfection strategy, there is no room for compromise. At the same time, your disinfection system must be able to achieve results at a cost you can afford.

The Halo Disinfection System® is validated to **kill 99.9999% of C. difficile spores** and ensures a uniform kill rate throughout every room. **The system combines this industry leading sporicidal efficacy with low cost of ownership and ease of use to deliver the greatest value in whole room disinfection on the market today.**



## Emergency Vehicles

*\*\* HaloMist™ has demonstrated effectiveness against viruses similar to the 2019-nCoV Coronavirus on hard, nonporous surfaces. Therefore, HaloMist™ can be used against the 2019-nCoV Coronavirus, now called SARS-CoV-2, when used in accordance with the directions for use against Feline calicivirus (as a surrogate for Norovirus) and Minute virus of Mice (MVM) on hard, nonporous surfaces. Refer to the CDC -or- OIE website at <https://www.cdc.gov/coronavirus/for additional information>.*

## Halosil's Five Star, Six Log Approach to Whole Room Disinfection.

- ★ **Efficacy:** Kills 99.9999% of resistant *C. difficile* spores, even in areas beyond human or UV reach.
- ★ **Affordability:** Reduces operating costs with a maintenance-free design and affordable disinfectants.
- ★ **Safety:** Safely decomposes into water, oxygen and non-toxic silver complexes.
- ★ **Versatility:** Flexibility to position dispensing nozzle assembly to maximize aerosol coverage within vehicles and spaces.
- ★ **Usability:** Empowers staff with intuitive operator controls and a simple, no rinse and no wipe formula.

## Halo Disinfection System®

# The Dry Fogging Technology Redefining Emergency Vehicle Disinfection.

### The Formula


Using hydrogen peroxide and ionic silver, Halosil's proprietary **HaloMist™** formula applies two mechanisms of killing action to destroy a broad spectrum of microorganisms and provide antimicrobial qualities for maximum impact.

### The Delivery

The **HaloFogger®** generates a turbulent aerosol that uses evaporation to ensure the uniform delivery of HaloMist™ throughout complex spaces—even in areas beyond the reach of humans and UV lights. HaloFogger's dry aerosol is safe for use around electronic equipment.

### The Results

The **Halo Disinfection System®** is validated to kill a diverse range of spores, bacteria, fungi and viruses. Specific kill claims can be found below.



## The HaloFogger EXT is the one complete system to handle all of your emergency vehicle (and facility!) disinfection needs.



The versatile **HaloFogger EXT** ensures complete surface disinfection by delivering the precise amount of HaloMist needed for treatment. The choice of a tripod or suction cup nozzle mount enables flexible nozzle placement to obtain maximum usability. Emergency vehicle interiors, locker rooms, equipment rooms, and social spaces can be disinfected with a single base HaloFogger unit. Additional nozzle assemblies can be purchased to treat multiple vehicles and spaces quickly and affordably.



**Halosil International**  
91 Lukens Drive  
New Castle, DE 19720 USA  
[www.halosil.com](http://www.halosil.com)



Phone:  
Website: